

NcMission™ hMSC Medium V3.0

Catalog # RP02010 1 Kit (525 mL)

Product Introduction

NcMission™ hMSC Medium V3.0 is a serum-free, animal component-free complete medium designed for the culture of primary human mesenchymal stem cells (hMSCs). hMSCs cultured in this medium exhibit stable proliferation, maintain normal expression of surface markers (CD73+/CD90+/CD105+, CD14-/CD34-/CD45-/CD79α-/HLA-DR-), and retain their trilineage differentiation potential (osteogenic, chondrogenic, and adipogenic differentiation).

Product Information

Table 1. NcMission™ hMSC Medium V3.0 Product Description

Product Information	Cat.No.	Amount	Storage
NcMission™ hMSC Medium V3.0 contains:	RP02010	1 Kit	2 °C~8 °C*
NcMission™ hMSC Medium V3.0 Basal Medium	RP02010-1	500 mL	2 °C~8 °C
NcMission™ hMSC Medium V3.0 Supplement (21×)	RP02010-2	25 mL	-20 °C or -80 °C

*Mix the basal medium and supplements to prepare the complete medium, which is stable at 2°C to 8°C for up to 2 weeks.

Reagents and Materials

Table 2. Reagents & Materials

Reagents & Materials	Brand (e.g.)	Cat.NO. (e.g.)
NcMission™ hMSC Medium V3.0	Shownin	RP02010
hMSC Cryopreservation Medium	Shownin	SN-06-1310
TrypLE Express Enzyme (1X), no phenol red	Thermo Sci.	12604013
T75/T175/T225 Culture Flasks	Thermo Sci.	156499/159910/159934
15 mL/50 mL Centrifuge Tube	Thermo Sci.	N/A
1.5/2 mL Cryovials	Thermo Sci.	N/A
10 µL/200 µL/1000 µL Pipette Tips	Rainin.	N/A

Preparation of Complete Medium

- 4.1 Thaw NcMission™ hMSC Medium V3.0 Supplement (21×) at 4 °C, and **do not thaw at 37°C**.
- 4.2 In a biosafety cabinet, use a sterile pipette to mix the following two components to prepare the complete culture medium.

NcMission™ hMSC Medium V3.0 Basal Medium: 500 mL

NcMission™ hMSC Medium V3.0 Supplement: 25 mL

- 4.3 The complete medium can be stored at 2 °C – 8 °C and should be used within 2 weeks.

Tips: The Supplement can be aliquoted based on actual requirements and stored frozen. For instance, it can be aliquoted into 5 vials, each containing 5 mL. Prior to use, thaw 5 mL of the supplement and combine it with 100 mL of basal medium to formulate the complete medium, which must be used within 2 weeks. The supplement should not be subjected to more than 2 freeze-thaw cycles.

Isolation and Culture of Primary hMSCs (Using Adipose Tissue as an Example)

- 5.1 Adipose Tissue Collection: Collect adipose tissue according to clinic protocols. Transport it at 4°C and process within 24 hours.
- 5.2 Under aseptic conditions, aspirate the lipoaspirate. Wash the tissue several times with DPBS (or physiological saline) to remove drugs used during liposuction and blood cells until the wash solution is no longer blood-colored. Using sterile ophthalmic scissors and forceps, clean the tissue and mince it into pieces of approximately 1-2 mm³.
- 5.3 Digest with 0.1% Type II collagenase at 37°C with agitation for 45-60 minutes. Centrifuge the digest at 800 × g for 10 minutes. The upper layer consists of undigested adipose tissue and oil. Carefully insert a pipette into the lower layer and aspirate the cell-containing digestate. Filter the digestate through a 70 µm cell strainer. Centrifuge the filtrate at 600 × g for 8 minutes. Discard the supernatant. Resuspend the cell pellet in a volume of DPBS (or physiological saline) equivalent to twice the pellet volume, and centrifuge at 600 × g for 5 minutes. Repeat this centrifugation wash step under the same conditions once more.

- 5.4 Add NcMission hMSC complete medium, adjust the cell density to 2×10^4 cells/mL, inoculate into T25 cm² culture flasks, and incubate in a 37°C, 5% CO₂ incubator.
- 5.5 Replace with fresh medium after 48 hours, and then change the medium every 3 days thereafter. Observe and record the growth status. When the cultured cells reach 90% confluence, perform subculture.
- 5.6 Cell digestion: Aspirate the culture supernatant and tissue pieces, add normal saline for washing once and then aspirate and discard. Add TrypLE pre-warmed to 37°C, with the amount of digestion solution. Digest at 37°C for 5-7 minutes (do not move midway), then add an equal volume of NcMission hMSC complete medium (or normal saline) to terminate digestion. Collect the cells and centrifuge (200 × g, 5 minutes).
- 5.7 Cell counting: Resuspend the cells in 5-10 mL of normal saline, filter through a 100 μm cell sieve once, and take a sample for counting. The cell viability should be ≥ 90%. Collect the cells by centrifugation (200 × g, 5 minutes).
- 5.8 Cell Seeding: Resuspend the cells in 5 mL of NcMission hMSC complete medium and seed at a density of 6000-8000 cells/cm². Add pre-warmed NcMission hMSC complete medium. Gently rock the flask horizontally three times and incubate in a 37°C, 5% CO₂ incubator. Shake again and proceed with culture. Culture for 3 days until 80-85% confluent before passaging.
- 5.9 Cell Cryopreservation: If cell cryopreservation is required, after centrifugation in step 5.6, add cryopreservation solution to resuspend the cells at a certain density (e.g., 2×10^6 cells per vial). Transfer to a programmable cooling box, freeze at -80°C overnight, and store in liquid nitrogen the next day.

Thawing hMSCs (Using a T75 Flask as an Example; the Procedure is Also Applicable to Other Culture Vessels)

- 6.1 Preheat a water bath to 37°C. Warm an appropriate amount of NcMission hMSC complete medium to room temperature.
- 6.2 Retrieve cryopreserved cells from dry ice and thaw in a 37°C water bath until only small ice crystals remain.
- 6.3 Transfer the cell suspension to a 15 mL centrifuge tube. Add 10 mL of pre-warmed NcMission hMSC

complete medium dropwise and mix gently. Centrifuge at $200 \times g$ for 5 minutes. Resuspend the cells in 5 mL of NcMission hMSC complete medium and count.

- 6.4 Seed the cells into the cell culture vessel at an appropriate seeding density (6000-8000 cells/cm²), and add an appropriate amount (refer to Table 3) of NcMission hMSC complete medium that has been equilibrated to room temperature. Gently rock the flask horizontally three times and incubate in a 37°C, 5% CO₂ incubator. Culture for 3 days until 80-85% confluent before passaging.

Table 3. Recommended Reagent Usage for hMSC Passaging and Culture

Culture Vessel	Growth Area	NcMission hMSC Complete Medium	TrypLE
6-Well Plate	9.6 cm ² per well	2 mL/well	1 mL/well
T75 Flask	75 cm ²	15 mL	4 mL
T175 Flask	175 cm ²	25 mL	8 mL
T225 Flask	225 cm ²	35 mL	10 mL

Passaging and Cryopreservation of hMSCs (Using T75 Flask as an Example; the Procedure is Also Applicable to Other Culture Vessels)

- 7.1 Selection of Passaging Timing: The growth rates of different hMSCs vary, and it is recommended to select the appropriate passaging time based on cell confluence. Passage cells when they reach 80-85% confluence.
- 7.2 Warm NcMission hMSC complete medium and digestion solution (Research-grade culture: Trypsin solution + Trypsin inhibitor; Clinical-grade culture: TrypLE) to room temperature.
- 7.3 Cell Detachment: Aspirate the medium, wash with DPBS (without Ca²⁺ /Mg²⁺), and add pre-warmed (37°C) TrypLE (referring to Table 3) . Incubate at 37°C for 5-7 minutes (do not move midway). Neutralize with an equal volume of NcMission hMSC complete medium. Centrifuge at $200 \times g$ for 5 minutes.
- 7.4 Cell Counting: Resuspend the cells in 5 mL of saline, filter through a 100 µm strainer, and count. Cell viability should be $\geq 90\%$. Centrifuge again at $200 \times g$ for 5 minutes.

7.5 Cell Seeding: Resuspend the cells in 5 mL of NcMission hMSC complete medium and seed at 6000-8000 cells/cm². Add pre-warmed fresh NcMission hMSC complete medium (referring to Table 3) . Gently rock the flask horizontally three times and incubate in a 37°C, 5% CO₂ incubator. Culture for 3 days until 80-85% confluent before passaging.

7.6 Cell Cryopreservation: If cell cryopreservation is required, after step 7.3, add cryopreservation solution to resuspend the cells at a certain density (e.g., 2×10^6 cells per vial). Transfer to a programmable cooling box, freeze at -80°C overnight, and store in liquid nitrogen the next day.

Adaptation of hMSCs from Other Culture Systems to NcMission™ hMSC Medium V3.0

When transitioning to NcMission™ hMSC Medium V3.0, recover or passage the cells in their original medium. On day 1, replace the medium with NcMission™ hMSC Medium V3.0. After one passage, the cells will adapt to the new system.