

ncMission™ M30 hMSC Medium (Phenol Red-free)

Product Manual

I. Product Introduction

ncMission™ M30 hMSC Medium (Phenol Red-Free) is a serum-free, animal-free complete medium developed by Anhui Shownin Biotech Co., Ltd. for the expansion of primary Human Mesenchymal Stem Cells (hMSCs). hMSCs cultured in this medium exhibit stable proliferation, maintain normal expression of cell surface markers (CD73+/CD90+/CD105+, CD14-/CD34-/CD45-/CD79α-/HLA-DR-), and retain their trilineage differentiation potential (osteogenic, chondrogenic, and adipogenic differentiation).

II. Product Information

Table 1: ncMission® M30 hMSC Amplification Medium (Phenol Red-free) Product Description

Product Information	Cat.No.	Amount	Storage
NcMission™ M30 hMSC Medium (Phenol Red-free) contains:	SN-02-0030	1 Kit	
NcMission™ M30 hMSC Medium Basal Medium (Phenol Red-free)	SN-00-0003	500 mL	2~8℃
NcMission™ M30 hMSC Medium serum-free Supplement	SN-02-0031	25 mL	-20℃ or -80℃

* Mix the base medium and supplements to prepare the complete medium. It can be stored at 2℃ to 8℃ and should be used within two weeks.

III. Reagents and Materials

Table 2: Reagents & Materials

Reagents & Materials	Brand (e.g.)	Cat.No (e.g.)
NcMission™ M30 hMSC Medium (Phenol Red-free)	Shownin	SN-02-0030
hMSC High-Efficiency Cryopreservation Medium	Shownin	SN-06-1310
hMSC Mild Digestion Solution	Shownin	SN-02-0040
T75/T175/T225 Culture Flasks	Thermo Sci.	156499 /159910/159934
15 mL/50 mL Centrifuge Tubes	Thermo Sci.	N/A
1.5/2 mL Cryovials	Thermo Sci.	N/A
10 µl/200 µl/1000 µl Pipette Tips	Rainin.	N/A
Freezing Container	Thermo Sci.	5100-0001

IV. Complete Medium Preparation

1. Thaw NcMission™ M30 hMSC Serum-Free Supplement (21×) at 4℃, and **do not thaw at 37℃**.
2. In a biosafety cabinet, mix the following components using sterile pipettes to prepare the complete medium:
NcMission® M30 hMSC Basal Medium: 500 mL
NcMission® M30 hMSC Serum-Free Supplement: 25 mL
3. The complete medium can be stored at 2-8℃ and should be used within 2 weeks.

TIPS: The serum-free supplement can be aliquoted and stored frozen according to usage. For example, aliquot 5 mL × 5 vials. Before use, thaw 5 mL of the supplement and mix with 100 mL of basal medium to prepare the complete medium, which should be used within 2 weeks. The supplement should not undergo more than two freeze-thaw cycles.

V. Isolation and Culture of Primary MSCs (Using Umbilical Cord Tissue Explant Method as an Example)

4. Umbilical Cord Collection: Collect the umbilical cord and place it in NcMission™ M30 hMSC Basal Medium. Transport at 4°C and process within 24 hours.
5. Material Preparation: Prepare freshly made NcMission™ M30 hMSC Medium, sterile culture dishes (6–10), medical-grade 75% ethanol, physiological saline, a tool kit (2 pairs of scissors, 2 pairs of forceps), and the umbilical cord in preservation solution. Transfer all materials to the biosafety cabinet.
6. Disinfection: Aspirate the preservation solution and submerge the umbilical cord in 75% ethanol for 2 minutes.
7. Washing: Transfer the umbilical cord to a sterile culture dish and wash 2–3 times with physiological saline to remove residual blood.
8. Cutting: Cut the umbilical cord into 2–3 cm segments and wash again with saline 2–3 times.
9. Separation: Cut along the vein to remove the venous wall, then remove the two arteries. Carefully separate Wharton's jelly, avoiding the epithelium.
10. Weighing: Transfer Wharton's jelly to a 50 mL centrifuge tube, add 3–5 drops of saline to keep it moist, and mince into 2–3 mm³ pieces. Weigh the tissue.
11. Seeding: Resuspend the minced tissue in NcMission™ M30 hMSC Medium and seed into culture flasks according to Table 3. Incubate at 37°C, 5% CO₂, and saturated humidity.
12. First Medium Change: On day 5, cells will begin to migrate out from the tissue explants. Tilt the culture flask upright at a 30-degree angle to allow the tissue explants to settle naturally in one corner of the flask. Aspirate the supernatant and slowly add fresh, pre-warmed NcMission™ M30 hMSC Medium. Gently mix and return the flask to the incubator for continued culture.
13. Second Medium Change: On days 9–10, the migrated cells will exhibit healthy morphology and begin to form multilayered growth. Tilt the culture flask upright at a 30-degree angle to allow the tissue explants to settle naturally in one corner of the flask. Aspirate the supernatant and slowly add fresh, pre-warmed NcMission™ M30 hMSC Medium. Gently mix and return the flask to the incubator for continued culture.
14. Passaging: Around day 12, cells can be passaged. Approximately 2–3×10⁶ cells/T75 flask (from 0.5 g Wharton's jelly) can be collected.
15. Cell Detachment: Aspirate the medium and tissue explants, wash once with saline, and add pre-warmed hMSC Mild Digestion Solution (refer to Table 4). Incubate at 37°C for 4–5 minutes, then add an equal volume of NcMission™ M30 hMSC Medium to neutralize the digestion. Centrifuge at 200×g for 5 minutes.
16. Cell Counting: Resuspend the cells in 5–10 mL of saline, filter through a 100 μm cell strainer, and count. Cell viability should be ≥90%. Centrifuge again at 200 × g for 5 minutes.
17. Cell Seeding: Resuspend the cells in 5 mL of NcMission™ M30 hMSC Medium and seed at a density of 5000–7000 cells/cm² (recommended 6000 cells/cm²). Add pre-warmed medium according to Table 4. Gently rock the flask horizontally three times, then incubate at 37°C, 5% CO₂, and saturated humidity. Rock the flask again three times and culture for 3 days. Passage when cells reach 80–85% confluence.
18. Cryopreservation: If cryopreservation is required, resuspend the cells in cryopreservation medium at a density of 2 × 10⁶ cells/mL after centrifugation. Transfer to a programmable freezing container, freeze at -80°C overnight, and then store in liquid nitrogen.

Table 3: Recommended Reagent Usage for Primary MSC Isolation Using Tissue Explant Method

Step	T75 Culture Flask	T175 Culture Flask	T225 Culture Flask
Wharton's Jelly	0.5 g	1 g	1.5 g

Weight			
Medium for Seeding	10 mL	15 mL	20 mL
1st Medium change (DAY5)	13 mL	20 mL	30 mL
2nd Medium change (DAY 9-10)	15 mL	25 mL	35 mL

VI. Passaging and Cryopreserving hMSCs (Using a T75 Flask as an Example; the Procedure is Also Applicable to Other Culture Vessels)

1. Preheat a water bath to 37°C. Allow ncMission® M30 MSC Expansion Medium to reach room temperature.
2. Retrieve the cryopreserved cells and transport them on dry ice to the cell culture room. Remove the cells from the dry ice and place them in a 37°C water bath. Gently agitate the vial until the ice crystals in the cell suspension are nearly completely dissolved (with only small, pea-sized ice crystals remaining). Remove the vial at this point.
3. Transfer the cell suspension to a 15 mL centrifuge tube and add 10 mL of room temperature NcMission™ M30 hMSC Medium dropwise. Centrifuge at 200 × g for 5 minutes, aspirate the supernatant, and resuspend in 5 mL of medium. Count the cells.
4. Seed the cells at 5000–7000 cells/cm² (recommended 6000 cells/cm²) and add pre-warmed medium according to Table 4. Rock the flask horizontally in a cross pattern three times, incubate at 37°C, 5% CO₂, and saturated humidity. Rock the flask again three times and culture for 3 days. Passage when cells reach 80–85% confluence.

Table 4: Recommended Reagent Usage for hMSC Passaging and Culture

Culture Vessel	Growth Area	NcMission™ Complete Medium	Trypsin/Trypsin Inhibitor
6-Well Plate	9.6 cm ² /hole	2 mL/well	1 mL/well
T75 Culture Flasks	75 cm ²	15 mL	4 mL
T175 Culture Flasks	175 cm ²	25 mL	8 mL
T225 Culture Flasks	225 cm ²	35 mL	10 mL

VII. Passaging & Cryopreserving hMSCs (Using a T75 Flask as an Example; the Procedure is Also Applicable to Other Culture Vessels)

5. **Passaging Timing:** The growth rate of different hMSC varies, and it is recommended to choose the exact timing of passaging based on the cell confluence, which can be done at around **80-85%** cell confluence.
6. **Preparation:** Allow NcMission™ M30 hMSC Medium and hMSC Mild Digestion Solution to reach room temperature.
7. **Cell Detachment:** Aspirate the medium, wash with DPBS (without Ca²⁺/Mg²⁺), and add pre-warmed digestion solution (refer to Table 4). Incubate at 37°C for 4–5 minutes, then neutralize with an equal volume of NcMission™ M30 hMSC Medium. Centrifuge at 200×g for 5 minutes.
8. **Cell Counting:** Resuspend the cells in 5 mL of saline, filter through a 100 μm cell strainer, and count. Cell viability should be ≥90%. Centrifuge again at 200×g for 5 minutes.
9. **Cell Seeding:** Resuspend the cells in 5 mL of NcMission™ M30 hMSC Medium and seed at 5000–7000 cells/cm² (recommended 6000 cells/cm²). Add pre-warmed medium according to Table 4. Rock the flask horizontally in a cross pattern three times, incubate at 37°C, 5% CO₂, and saturated humidity. Rock the flask again three times and culture for 3 days. Passage when cells reach 80–85% confluence.
10. **Cryopreservation:** Resuspend the cells in cryopreservation medium at 2×10⁶ cells/mL after centrifugation. Transfer to a programmable freezing container, freeze at -80°C overnight, and then store in liquid nitrogen.

VIII. Adaptation of hMSCs from Other Culture Systems to NcMission™ M30 hMSC Medium

When transitioning hMSCs from other culture systems to NcMission™ M30 hMSC Medium, recover or passage the cells in the original medium, then switch to NcMission™ M30 hMSC Medium on day 1. The cells should adapt to the new medium after one passage.