

NcTarget™ hPSC Medium

Product Manual

Catalog#RP01020 500mL

I. Product Introduction

NcTarget™ hPSC Medium is a serum-free, chemically defined medium designed for use under feeder-free condition with human pluripotent stem cells (hPSCs), including human embryonic stem cells (hESCs) and human induced pluripotent stem cells (hiPSCs). hPSCs proliferate rapidly in NcTarget™ hPSC Medium, while differentiated cells grow more slowly, enabling selective expansion and acquisition of high purity hPSC.

II. Product Information

Table 1: NcTarget™ hPSC Medium Product Description

Product Information	Cat.NO.	Amount	Storage
NcTarget™ hPSC Medium Contains:	RP01020	1 Kit	2 °C~8 °C*
NcTarget™ hPSC Basal Medium	RP01020-1	400 mL	2 °C~8 °C
NcTarget™ hPSC Supplement-A	RP01020-A	20 mL	-20 °C or -80 °C
NcTarget™ hPSC Supplement-B	RP01020-B	80 mL	-20 °C or -80 °C

***Mix the basal medium and supplements to prepare the complete medium, which can be stored at 2 °C to 8 °C and used within 2 weeks.**

III. Reagents and Materials

Table 2: Reagents & Materials

Reagents & Materials	Brand (e.g.)	Cat.NO. (e.g.)
NcTarget™ hPSC Medium	Shownin	RP01020
hPSC Cryopreservation Medium	Shownin	RP01003
hPSC Dissociation Buffer	Shownin	RP01007
Blebbistatin (10 mM)	Shownin	RP01008
DMEM/F12 Medium	Thermo Sci.	11330
6-well Plate	Thermo Sci.	140685
1 mL/5 mL/10 mL/25 mL Pipettes	Thermo Sci.	N/A
15 mL/50 mL Centrifuge Tubes	Thermo Sci.	N/A
1.5/2 mL Cryovials	Thermo Sci.	N/A
Freezing Container	Thermo Sci.	5100-0001

IV. Reagent Preparation

(i) Preparation of Complete NcTarget™ hPSC Medium (500 mL)

1. Thaw **NcTarget™ Supplement A/B** at 4 °C; **do not thaw at 37 °C**.
2. In a biosafety cabinet, mix the following components using sterile pipettes to prepare the complete medium.
The complete medium can be stored at 4 °C and used within 2 weeks.

Table 3: NcTarget™ hPSC Medium-Preparation Guide*

Component	Final Volume		
	500 mL	100 mL	50 mL
NcTarget™ Basal Medium	400 mL	80 mL	40 mL
NcTarget™ Supplement-A	20 mL	4 mL	2 mL
NcTarget™ Supplement-B	80 mL	16 mL	8 mL

*Aliquot and freeze NcTarget™ Supplement A/B based on actual usage. Do not exceed 2 freeze-thaw cycles for NcTarget™ Supplement A/B.

(ii) Matrigel® Coating of Culture Plates (Example: 6-well plates coating with Corning® Matrigel®)

A. Aliquoting Matrigel®

1. Determine the concentration of Matrigel® based on the batch number. Calculate the aliquot volume and quantity based on the recommended coating concentration (0.013 mg/cm² for hPSC culture).
For example, if Matrigel® concentration is 11.3 mg/mL, aliquot 3 mg per tube (sufficient for 4 six-well plates).
2. Prepare 38 sterile 1.5 mL EP tubes, label properly. Pre-cool tubes and pipette tips at -20 °C for 1 hour.

TIPS: For hESC-Qualified Matrigel® (Cat. No. 354277), follow the Dilution Factor provided (e.g., 238 µL for 4 plates).

$$\text{Number of Aliquots} = 5 \text{ mL} / 238 \text{ µL} = 21.01.$$

3. Thaw Matrigel® at 4 °C overnight. Ensure Matrigel® is fully thawed before aliquoting.
TIPS: Matrigel® remains liquid only at 4 °C. Temperature fluctuations may cause it to solidify.
4. Place thawed Matrigel®, pre-cooled EP tubes, and pipette tips on ice.
5. Mix Matrigel® and aliquot into pre-cooled EP tubes. Replace pipette tips if clogged to ensure accurate volumes.
6. Store aliquots at -20 °C.

B. Plating

1. Prepare 36 mL of chilled DMEM/F12 in a 50 mL centrifuge tube.
2. Pre-cool 1000 µL pipette tips at -20 °C for 1 hour. Thaw one Matrigel® aliquot (3 mg) at 4 °C.
3. Place thawed Matrigel® and pre-cooled pipette tips on ice.
4. Add 1 mL of cold DMEM/F12 to the thawed Matrigel® (3 mg). Gently pipette to mix thoroughly.
5. Transfer the Matrigel® mixture to the remaining DMEM/F12 in the centrifuge tube. Mix well using a 10 mL pipette.
6. Add 1.5 mL/well to 4 six-well plates. Gently rock to ensure even coating.

- Let plates sit at room temperature for 1 hour or at 4 °C overnight for immediate use. Or it can be stored at 4 °C for use within 2 weeks.

V. Thawing and Recovery of hPSC (Using 6-well plate operation as an example, the procedure is also applicable to other vessels)

- Preheat the water bath to 37 °C.
- Place a Matrigel®-coated 6-well plate in the biosafety cabinet for 1 hour to warm to room temperature (15–30 °C).
- Prepare 4 mL of **NcTarget™ Complete Medium** with 1 µL Blebbistatin (10 mM) at a 1:4000 ratio. Warm to room temperature (15–30 °C).

TIPS: Do not pre-warm the medium in a 37 °C water bath.

- Thaw a frozen cell vial in the 37 °C water bath, gently shaking until ice crystals nearly disappear (within 1 minute).
- Wipe the vial with 75% ethanol and transfer to the biosafety cabinet. Transfer the cell suspension to a 15 mL centrifuge tube. Slowly add 10 mL of DMEM/F12, mix gently, and centrifuge at 160 × g for 5 minutes.
- Aspirate the supernatant. Resuspend the cell pellet in 4 mL of pre-warmed Blebbistatin + **NcTarget™ Complete Medium**. Avoid excessive pipetting.
- Aspirate the Matrigel® coating solution from two wells of the 6-well plate. Seed the cell suspension at 2 mL/well.
- Gently rock the plate horizontally in a cross pattern 3 times. Place in a 37 °C, 5% CO₂ incubator with saturated humidity. Rock the plate 3 more times and maintain it under standard conditions overnight.
- Replace the medium with fresh **NcTarget™ Complete Medium** after 18–24 hours. Change the medium daily thereafter.

Table 3: Recommended Usage of Reagents

Culture Vessel	Growth Area	DPBS (mL)	hPSC Dissociation Buffer	hPSC medium*
6-well Plate	9.6 cm ² /well	2 mL/well	2 mL/well	2 mL/well
12-well Plate	4.5 cm ² /well	1 mL/well	1 mL/well	1 mL/well
24-well Plate	2 cm ² /well	0.5 mL/well	0.5 mL/well	0.5 mL/well
35mm Culture Dish	8 cm ²	2 mL	2 mL	2 mL

* We recommend adding 50% more medium when confluence exceeds 50%. For example, add 3 mL/well to a 6-well plate during medium change.

VI. Passaging of hPSC

(Using 6-Well Plates and EDTA Dissociation as an Example; the Procedure is Also Applicable to Other Vessels)




Figure 1: Morphology of hiPSCs under NcTarget™ hPSC medium use conditions on a Matrigel® Plate.

(A, B, C, D) Morphology on days 1, 2, 3, and 4, respectively. Scale bar: 200 μm.

1. Determining the Timing for Passaging:

Passage iPSCs when either of the following conditions are met:

- 1.1 When cell confluence reaches approximately 85% (Figure 1-C, D), normally 4 days after seeding. Even if the clones are small and confluence is insufficient, continuous culture should not exceed 5 days.
- 1.2 When cell confluence is low but the colonies are overly large, leading to poor growth in the central region of colonies.

2. Passaging ratio:

The passaging ratio can range from 1:5 to 1:20 depending on cell growth status and experimental requirements. If cells are healthy, with 85% confluence and uniform clone size (Figure 1-C, D), a 1:10 ratio is recommended. Reduce the ratio if cell density is low; increase the ratio if density is high.

3. Preparation of Vessels:

Place Vitronectin-coated 6-well plates in the biosafety cabinet for approximately 1 hour to equilibrate to room temperature (~25 °C).

4. Prepare 2 mL/well of **NcTarget™ complete medium** according to the number of wells for the next step. Add Blebbistatin (10 mM) at a 1:4000 ratio and equilibrate to room temperature (~25 °C).

Tip: Add 0.5 µL of Blebbistatin (10 mM) to 2 mL of NcTarget™ complete medium.

5. Aspirate the medium from the iPSC wells and add 2 mL/well of DPBS (without calcium and magnesium). Gently swirl and aspirate.
6. Add 2 mL/well of **hPSC Dissociation Buffer** to fully cover the well surface.
7. Maintain in a 37 °C incubator for 7–8 minutes.

Tips: (1) Observe cell morphology under a microscope after 7–8 minutes. Stop dissociation when most cells appear bright and round but have not yet detached or floated (Figure 2C). If most cells remain unchanged, extend the dissociation time (Figures 2A & 2B).

(2) Place the 6-well plate directly on the metal shelf in the incubator for even heating. Do not stack plates.




Figure 2: (A) Digestion for 4 min; (B) Digestion for 6 min; (C) Digestion for 8 min. Scale bar: 200 µm.

8. Bring the plate to the biosafety cabinet without shaking. Tilt the plate and aspirate the hPSC Dissociation Buffer (0.5 mM EDTA).
9. Immediately add 2 mL/well of pre-warmed Blebbistatin-supplemented **NcTarget™ complete medium**. Gently rock the plate in a horizontal cross pattern to detach cells.

Tips: (1) Gently pipette 1–2 times after adding the medium. Avoid excessive pipetting (>2 times).

(2) Avoid scraping cells. It is normal for 10–15% of cells to remain attached. If a large number of cells remain attached, extend the dissociation time.

(3) Process no more than one 6-well plate at a time. Quickly aspirate the medium after adding NcTarget™ medium, as EDTA dissociation is rapidly neutralized and the cells re-attach quickly in medium. hPSCs should not remain in hPSC Dissociation Buffer (0.5 mM EDTA) for >15 minutes.

10. Seeding:

- 10.1 Aspirate the Vitronectin solution from the 6-well plate and add 2 mL/well of pre-warmed Blebbistatin + **NcTarget™ Complete Medium**.

- 10.2 Gently mix the cell suspension obtained in step 9 and evenly distribute it into the wells according to the predetermined passage ratio.

Tips: Alternatively, calculate the total cell number required per plate, transfer the suspension accordingly to a 15 mL centrifuge tube, and bring up the volume to 12 mL with pre-warmed Blebbistatin-supplemented NcTarget™ complete medium. Distribute evenly into Matrigel®-coated 6-well plates after aspirating the coating solution.

11. Gently rock the 6-well plate horizontally in a cross pattern 3 times. Place it in a 37 °C, 5% CO₂, humidified incubator. Rock the plate again 3 times and maintain the plate under standard conditions overnight.
12. After 18–24 hours, replace with fresh **NcTarget™ Complete Medium**. Change the medium daily thereafter, followed with passaging or cryopreservation 4–5 days later (Figure 1).

VII. Cryopreservation of hPSC

1. When the cell confluence reaches about 85% (Figures 1), cells can be harvest for cryopreservation. Typically, 2×10^6 – 4×10^6 viable cells per well of a 6-well plate can be collected and frozen in one cryovial.
2. Prepare the appropriate number of 1.5/2 mL cryovials and label properly.
3. Take out hPSC High-Efficiency Cryopreservation Medium from the 4 °C refrigerator and equilibrate to room temperature. Mix thoroughly before use.

Tips: DMSO in the cryopreservation medium tends to settle at the bottom. Inadequate mixing may result in insufficient DMSO concentration initially and excessive concentration later, leading to instability in cryopreserved cells.

4. Aspirate the medium from the plate and add 2 mL/well of DPBS (without calcium and magnesium). Gently rock swirl several times and aspirate.
5. Add 2 mL/well of hPSC dissociation working solution. Place the plate in a 37 °C incubator for 7–8 minutes (refer to Section VI, Step 7).
6. After digestion, gently take out the culture plate and aspirate the EDTA buffer.
7. Mix the pre-warmed hPSC cryopreservation medium thoroughly. Add 1 mL of cryopreservation medium to each well. Gently pipette to resuspend the cells, then rock the plate in a horizontal cross pattern 3 times. Transfer the cell suspension into the labeled 1.5/2 mL cryovials.
8. Place the vials in a cell freezing container and store at -80 °C overnight. Transfer the vials to a liquid nitrogen tank for long-term storage the next day. Alternatively, use a programmable controlled rate freezer to cool the cells below -80 °C before transferring to liquid nitrogen storage.

VIII. Adaptation of hPSC from Other Medium Systems to NcTarget™ hPSC Medium

hPSCs under other feeder-free medium system can be transitioned to **NcTarget™ complete medium** when the cells are in good condition. Initially, replace the medium with a 1:1 mixture of the original medium and **NcTarget™ complete medium** 24 hours after passaging with the original medium. Use the 1:1 mixture for the next two passages before completely switch to **NcTarget™ medium**. The cells need 2~3 passages in NcTarget™ medium to fully adapt to the new medium system. Cells can be cryopreserved after fully adaptation.

IX. Problems and Solutions

➤ Differentiation Observed in hiPSC handling
<ul style="list-style-type: none"> • Ensure NcTarget™ complete medium is stored at 4 °C and used within 2 weeks. Only pre-warm the amount of medium required for the current experiment to minimize temperature fluctuations and prevent degradation of ingredients in the medium. • If hPSC colonies exhibit good overall morphology with sporadic differentiated cells (<1%) at the edges, these cells can be removed during EDTA passaging. • Ensure that the size of the hPSC cell clusters during passaging is uniform, with the size of about 20 cells being ideal. For larger clusters, gently pipette no more than 3 times using a 5 mL pipette, applying light and even pressure to avoid cell damage or differentiation. • Avoid removing cells from the incubator for more than 15 minutes during observation. • If hPSC clones appear loose internally with irregular edges and differentiation exceeds 20%, discard.
➤ Can Dispase or Collagenase Be Used for hPSC Passaging?
<ul style="list-style-type: none"> • Dispase or collagenase can be used, but digestion efficiency may be suboptimal, affecting post-passaging cell viability and potentially accumulating differentiated cells.. • For hPSCs cultured in the NcTarget™ system, non-enzymatic, gentle dissociation methods are recommended for passaging. • If single-cell dissociation is required for experiments, use Accutase for 5–10 minutes.
➤ hPSCs Fail to Adhere or Exhibit Low Adhesion After Passaging
<ul style="list-style-type: none"> • Avoid excessive passaging ratios (>1:20). • EDTA dissociation time should not be too long. While some cell lines may require extended dissociation beyond 8 minutes, do not exceed 15 minutes. • Avoid excessive pipetting (>3 times) to prevent clusters disruption or cell damage. • Ensure culture plates are coated with Vitronectin, Matrigel® • , or other matrices suitable for pluripotent stem cell growth. • ROCK inhibitor is necessary in the medium after passaging.
➤ Cells detach after medium change
<ul style="list-style-type: none"> • Perform the first medium change 18–24 hours after seeding to ensure proper cell adhesion. • Handle medium changes gently to avoid dislodging cell clusters from the matrix. • If cell seeding density is very low (e.g., for cell cloning experiments), avoid medium changes for 2–3 days, ensuring the medium contains ROCK inhibitor.
➤ Uneven distribution of hiPSC clusters in wells
<ul style="list-style-type: none"> • Ensure the coating matrix is evenly distributed across the bottom of the culture vessel. • During passaging, ensure cells are evenly dispersed. After rocking the plate in a horizontal cross pattern, avoid moving the plate to prevent cell aggregation in the center. • After placing the plate in the incubator, rock it again in a horizontal cross pattern to ensure even distribution.